AB-NP044-1 SHIP1-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human SHIP1 (INPP5D)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Phosphatidylinositol-3,4,5-trisphosphate 5-phosphatase 1
Alias:	Inositol polyphosphate-5-phosphatase of 145 kDa; Inositol polyphosphate-5-phosphatase, 145kDa; INPP5D; P150Ship; Phosphatidylinositol-3,4,5-trisphosphate 5-phosphatase 1; SH2-containing inositol 5'-phosphatase 1; SHIP1; SHIP-1; SIP-145
UniProt ID:	Q92835
Sequence Predicted Mass (KDa):	133.292 (1.189 AA; Q92835-1); 133.193 (1.188 AA; Q92835-2); 109.294 (976 AA; Q92835-3)
Observed SDS-PAGE Mass (KDa):	132-145

Immunogen

Antibody Immunogen Source:	Human SHIP1 (INPP5D) sequence peptide Cat. No.: PE-01BBP99
Antibody Immunogen Sequence:	CSEKEWLEILKHSLQE
Location in Target:	Corresponds to amino acid residues S457 to E471
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NP044-1 SHIP1-1 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

I I I I I I I I I I I I I I I I I I I	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	1 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Very high selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.