

# AB-PK511

## ACK1-pY284 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase ACK1 (TNK2)



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Activated CDC42 kinase 1
<b>Alias:</b>	ACK; Activated Cdc42-associated kinase 1; Activated p21cdc42Hs kinase; Activated p21cdc42Hs kinase 1; Kinase ACK1; Non-receptor protein tyrosine kinase Ack; p21cdc42Hs; TNK2; Tyrosine kinase, non-receptor, 2; FLJ44758; FLJ45547; ENSG00000061938
<b>UniProt ID:</b>	Q07912
<b>Sequence Predicted Mass (KDa):</b>	119.349 (1086 AA; Q07912-3); 114.569 (1038 AA; Q07912); 60.062 (528 AA; Q07912-2)
<b>Observed SDS-PAGE Mass (KDa):</b>	120-130

### Immunogen

<b>Antibody Immunogen Source:</b>	Human ACK1 (TNK2) sequence peptide Cat. No.: PE-04AIA80
<b>Antibody Immunogen Sequence:</b>	DDH(pY)VMQ(bA)C (bA) = beta-alanine
<b>Location in Target:</b>	Corresponds to amino acid residues D281 to Q287; In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.
<b>Peptide Type:</b>	For phosphosite-specific recognition of target.
<b>Target Phosphosite:</b>	Tyr-284

### Production

<b>Antibody Host Species:</b>	Rabbit
<b>Antibody Type:</b>	Polyclonal
<b>Antibody Ig Isotype Clone Lot:</b>	Immunoglobulin G
<b>Production Method:</b>	The immunizing peptide was produced by solid phase synthesis on a multipепeptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
<b>Antibody Amount:</b>	25 µg
<b>Antibody Concentration:</b>	1 mg/ml
<b>Lot Number:</b>	150203
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerazol
<b>Storage Conditions and Stability:</b>	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-PK511

## ACK1-pY284 Antibody



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Applications

<b>Product Use:</b>	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 µg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals
<b>Antibody Positive Controls:</b>	Very strong immunoreactivity with immunogen peptide on dot blots.
<b>Detection by Immunoblotting in Cell/Tissue Lysates:</b>	Very weak immunoreactivity of a target-sized protein by Western blotting in FBS-treated A431 cells.
<b>Overall Antibody Specificity:</b>	High selectivity
<b>Antibody Cross Reactivities:</b>	Cross-reactivity with CDK1 and ERK1 in seastar oocytes, and with a 27 KDa protein in A431 and HEK-293 cells.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINEXUS(546-3987)