## **AB-PK517**

Akt1-pY326 Antibody Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase Akt1 (PKBa)



Email: info@kinexus.ca Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	RAC-alpha serine/threonine-protein kinase
Alias:	AKT; Akt1; AKT1 kinase; C-AKT; Kinase Akt1; PKB; PKB-alpha; PRKBA; V-akt murine thymoma viral oncogene 1; Protein kinase B; RAC; RAC-alpha serine,threonine kinase; RAC-alpha serine/threonine kinase; RAC-PK-alpha; PKBA; hCG_96740; AKT1_NEW; AKT; MGC99656; CCDS9994; B0LPE5; B2RAM5; B3KVH4; B3KXD7; B7Z5R1; Q9BV07; ENSG00000142208
UniProt ID:	P31749
Sequence Predicted Mass (KDa):	55.686 (480 AA; P31749); 48.347 (418 AA; P31749-2)
Observed SDS-PAGE Mass (KDa):	57-62
Immunogen	
Antibody Immunogen Source:	Human Akt1 (PKBa) sequence peptide Cat. No.: PE-04AID95
Antibody Immunogen Sequence:	DND(pY)GRA(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues D323 to A329; In the protein kinase catalytic domain.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-326
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 µg
Antibody Concentration:	0.75 mg/ml
Lot Number:	141024
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-PK517 Akt1-pY326 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other vertebrates; Phosphosite is highly conserved in diverse species
Antibody Positive Controls:	Medium immunoreactivity with immunogen peptide on dot blots.
<b>Overall Antibody Specificity:</b>	Very high selectivity
Antibody Cross Reactivities:	Weak immunoreactivity on protein dot blots with human recombinant Akt2. No cross-reactive proteins detected in HepG2cells, MCF7 cells and seastar oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.