AB-PK547 BRK-pS446+pY447 Antibody

Torget Dretain

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase BRK



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca	
Phone: 604-323-2547	

Target Protein	
Name Long:	Protein-tyrosine kinase 6
Alias:	Brk6; PTK6; Tyrosine-protein kinase 6; FLJ42088; CCDS13524.1; Q13882; ENSG00000101213
UniProt ID:	Q13882
Sequence Predicted Mass (KDa):	51.834 (451 AA; Q13882); 14.465 (134 AA; Q13882-2)
Observed SDS-PAGE Mass (KDa):	50-55
Immunogen	
Antibody Immunogen Source:	Human BRK sequence peptide Cat. No.: PE-04AHP00
Antibody Immunogen Sequence:	SFT(pS)(pY)ENP(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues S443 to P450; In the C-terminus of the kinase
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Ser-446+Tyr-447
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	150305
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK547 BRK-pS446+pY447 Antibody

12 12



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium immunoreactivity with immunogen peptide on dot blots.
Detection by Immunoblotting in Cell/Tissue Lysates:	Weak immunoreactivity of a target-sized protein by Western blotting in A549 cells.
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No significant cross-reactivities detected in A549 except EGF increase in 75 KDa protein. Five higher molecular mass cross-reactive proteins (~75 to 250 KDa) in sea star oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)