# AB-PK570 CDK5-pY15 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase CDK5



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:	Cyclin-dependent protein-serine kinase 5
Alias:	Cell division protein kinase 5; CRK6; PSSALRE; TAU protein kinase II; TPKII; ENSG00000164885
UniProt ID:	Q00535
Sequence Predicted Mass (KDa):	33.304 (292 AA; Q00535); 29.544 (260 AA; Q00535-2)
Observed SDS-PAGE Mass (KDa):	29-35

### Immunogen

Antibody Immunogen Source:	Human CDK5 sequence peptide Cat. No.: PE-04AJO95
Antibody Immunogen Sequence:	EGT(pY)GTV(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues E12 to V18; In the protein kinase catalytic domain in subdomain I.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-15

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	0.5 mg/ml
Lot Number:	141204
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-PK570 CDK5-pY15 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other vertebrates; Phosphosite is highly conserved in diverse species
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in A431 cells and sea star oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.