## AB-PK590 DAPK-pS269 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase DLK (ZPK)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	Mitogen-activated protein kinase kinase kinase 12
Alias:	Dual leucine zipper bearing kinase; Dual leucine zipper kinase DLK; EC 2.7.11.25; Kinase DLK; Leucine-zipper protein kinase; M3K12; MAP3K12; MAPK-upstream kinase; MEKK12; MUK; ZPK; ZPKP1
UniProt ID:	Q12852
Sequence Predicted Mass (KDa):	96.322 (892 AA; Q12852-2); 93.219 (859 AA; Q12852)
Observed SDS-PAGE Mass (KDa):	
Immunogen	
Antibody Immunogen Source:	Human DLK (ZPK) sequence peptide Cat. No.: PE-04ALS85
Antibody Immunogen Sequence:	TKM(pS)FAG(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues T266 to G272; In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Ser-269
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	141024
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-PK590 DAPK-pS269 Antibody

12 12



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Very strong immunoreactivity with immunogen peptide on dot blots.
Detection by Immunoblotting in Cell/Tissue Lysates:	Very strong immunoreactivity of a target-sized protein by Western blotting in A431 cells.
<b>Overall Antibody Specificity:</b>	High selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in HepG2 cells and sea star oocytes, except for a ~76 KDa cross-reactive protein and phospho-CDK1 in germinal vesicle-positive, maturing sea star oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)