# **AB-PK626** ERK5-pY221 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase ERK5

(MAPK7)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:	Extracellular regulated protein-serine kinase 5 (Big MAP kinase 1 (BMK1))
Alias:	BMK; BMK1; BMK1 kinase; ERK4; ERK-5; Extracellular-signal-regulated kinase 5; Kinase ERK5; MAPK7; PRKM7; ENSG00000166484
UniProt ID:	Q13164
Sequence Predicted Mass (KDa):	88.386 (816 AA; Q13164); 73.218 (677 AA; Q13164-2); 59.328 (533 AA; Q13164-3); 50.152 (451 AA; Q13164-4)
Observed SDS-PAGE Mass (KDa):	88-115

### Immunogen

Antibody Immunogen Source:	Human ERK5 (MAPK7) sequence peptide Cat. No.: PE-04AEI99
Antibody Immunogen Sequence:	MTE(pY)VAT(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues M218 to T224; In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-221

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	141003
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-PK626 ERK5-pY221 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium immunoreactivity with immunogen peptide on dot blots.
Detection by Immunoblotting in Cell/Tissue Lysates:	Medium immunoreactivity of a target-sized protein by Western blotting in A431 cells.
Overall Antibody Specificity:	Medium-High selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in T98G cells, but manycross-reactive proteins detected in A431 cells, including phospho-ERK and phospho-p38 MAP kinases.

This product is for in vitro research use only and is not intended for use in humans or animals.