## AB-PK630 FER-pY402 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase Fer (TYK3)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	Fer (fps/fes related) protein-tyrosine kinase
Alias:	AV082135; C-FER; FER; Fer (fps/fes related) tyrosine kinase; FERT2; Kinase Fer; P94-FER; Phosphoprotein NCP94; TYK3; PPP1R74; CCDS4098.1; ENSG00000151422; B4DDX7
UniProt ID:	P16591
Sequence Predicted Mass (KDa):	94.638 (822 AA; P16591); 74.266 (647 AA; P16591-2); 51.643 (453 AA; P16591-3)
Observed SDS-PAGE Mass (KDa):	90-100
Immunogen	
Antibody Immunogen Source:	Human Fer (TYK3) sequence peptide Cat. No.: PE-04AMG99
Antibody Immunogen Sequence:	VVN(pY)EED(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues V399 to D405; In the region preceding the SH2 domain. One of the two major in vivo phosphorylations sites in Fer.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-402
Production Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The
	sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification
Antibody Amount: Antibody Concentration:	immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
	immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose. 25 μg
Antibody Concentration:	<ul> <li>immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.</li> <li>25 μg</li> <li>1 mg/ml</li> </ul>

## AB-PK630 FER-pY402 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Very strong immunoreactivity with immunogen peptide on dot blots.
<b>Overall Antibody Specificity:</b>	High-very high selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in sea star oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)