AB-PK722 Mos-pY263 Antibody

Target Protein

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase Mos



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	Moloney sarcoma oncogene-encoded protein-serine kinase
Alias:	ADRB2; C-mos; MGC119963; MSV; MGC119962; CCDS6164.1; ENSG00000172680
UniProt ID:	P00540
Sequence Predicted Mass (KDa):	37.820 (346 AA; P00540)
Observed SDS-PAGE Mass (KDa):	40-45
Immunogen	
Antibody Immunogen Source:	Human Mos sequence peptide Cat. No.: PE-04AOC99
Antibody Immunogen Sequence:	ADI(pY)SFA(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues A260 to A266; In the protein kinase catalytic domain near subdomain IX.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-263
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification
Antibody Amount:	over phosphotyrosine-agarose. 25 μg
Antibody Concentration:	0.5 / 1.0 mg/ml
Lot Number:	150305 / 160121
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK722 Mos-pY263 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium immunoreactivity with immunogen peptide on dot blots.
Overall Antibody Specificity:	High-very high selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in phenylarsine oxide (PAO)+vanadate-treated HeLa (homogenized in SDS-PAGE sample buffer), a strong ~160K cross-reactive protein in EGF-treated A431 cells (homogenized in SDS-PAGE sample buffer); In insulin-treated MCF7 cells, strong background >75 KDa was observed (homogenized in SDS-PAGE sample buffer).

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINEXUS(546-3987)