# AB-PK723 MSK1-pS212 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase MSK1 (RPS6KA5)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:	Mitogen & stress-activated protein-serine kinase 1; Ribosomal protein S6 kinase alpha 5
Alias:	90 kDa ribosomal protein S6 kinase, , polypeptide 5; 90 kDa ribosomal protein S6 kinase, KS6A5; Kinase MSK1; KS6A5; MGC1911; MSPK1; Nuclear mitogenand stress-activated protein kinase-1; Polypeptide V; Ribosomal protein S6 kinase; RLPK; RPS6KA5; MGC1911; MSPK1; RLPK; CCDS9893.1; ENSG00000100784
UniProt ID:	O75582
Sequence Predicted Mass (KDa):	89.865 (802 AA; O75582); 81.780 (723 AA; O75582-3); 61.769 (549 AA; O75582-2)
Observed SDS-PAGE Mass (KDa):	90-95

### Immunogen

Antibody Immunogen Source:	Human MSK1 (RPS6KA5) sequence peptide Cat. No.: PE-04AKQ95
Antibody Immunogen Sequence:	RAY(pS)FCG(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues R209 to G215; In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Ser-212

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	141202
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-PK723 MSK1-pS212 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium immunoreactivity with immunogen peptide on dot blots.
Overall Antibody Specificity:	Very high selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.