AB-PK749 PAK1-pT423 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase PAK1 (PAKa)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

	1	_	
l or	$\alpha \alpha t$	Uro	toin
		Pro ^o	
	gui		

Name Long:	p21-activated kinase 1 alpha; Protein-serine/threonine kinase PAK 1	
Alias:	ADRB2; Alpha-PAK; CDC42,RAC effector kinase PAK-A; CDC42/RAC effector kinase PAK-A; Kinase PAK1; MUK2; P21 protein (Cdc42/Rac)-activated kinase 1; P65-PAK; STE20, yeast; P68-PAK; PAK 1; PAK-1; PAKA; Protein kinase MUK2; PAKalpha; MGC130000; MGC130001; CCDS8250.1; Q13153; ENSG00000149269	
UniProt ID:	Q13153	
Sequence Predicted Mass (KDa):	61.632 (553 AA; Q13153-2); 60.647 (545 AA; Q13153)	
Observed SDS-PAGE Mass (KDa):	60-67	

Immunogen

Antibody Immunogen Source:	Human PAK1 (PAKa) sequence peptide Cat. No.: PE-04AOH95	
Antibody Immunogen Sequence:	KRS(pT)MVG(bA)C (bA) = beta-alanine	
Location in Target:	Corresponds to amino acid residues K420 to G426; In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.	
Peptide Type:	For phosphosite-specific recognition of target.	
Target Phosphosite:	Thr-423	

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	150305
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK749 PAK1-pT423 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays	
Antibody Dilution Recommended:	2 μg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals	
Antibody Positive Controls:	Very strong immunoreactivity with immunogen peptide on dot blots.	

This product is for in vitro research use only and is not intended for use in humans or animals.