AB-PK757 PDGFRa-pS847+pY849 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase PDGFRa



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Platelet-derived growth factor receptor kinase alpha
Alias:	Al115593; Alpha platelet-derived growth factor receptor precursor; CD140a; Kinase PDGFR-alpha; MGC74795; PDGFR2; PDGFRa; PDGFR-alpha; PDGF-R-alpha; PGFRA; Platelet-derived growth factor, alpha-receptor; Rhe-PDGFRA; RHEPDGFRA; CD140A; MGC74795; CCDS3495.1; ENSG00000134853
UniProt ID:	P16234
Sequence Predicted Mass (KDa):	122.670 (1089 AA; P16234); 82.809 (743 AA; P16234-3); 24.023 (218 AA; P16234-2)
Observed SDS-PAGE Mass (KDa):	115-135

Immunogen

Human PDGFRa sequence peptide Cat. No.: PE-04ARN60
D(pS)N(pY)VSKG(bA)C (bA) = beta-alanine
Corresponds to amino acid residues D846 to G853; In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
For phosphosite-specific recognition of target.
Ser-847+Tyr-849

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	141003
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK757 PDGFRa-pS847+pY849 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity with immunogen peptide on dot blots. Very strong immunoreactivity with recombinant human PDGFRA on protein dot blots.
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	Strong immunoreactivity on protein dot blots with recombinant human PDGFRB, FLT3 and Kit, and weak immunoreactivity with CSF1R.

This product is for in vitro research use only and is not intended for use in humans or animals.