AB-PK779 Plk1-pY217 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase Plk1

(PLK)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Polo-like protein-serine kinase 1
Alias:	ADRB2; Kinase PLK1; PLK; PLK-1; Polo-like kinase 1; Serine- threonine protein kinase 13; STPK13; CCDS10616; ENSG00000166851
UniProt ID:	P53350
Sequence Predicted Mass (KDa):	68.255 (603 AA; P53350)
Observed SDS-PAGE Mass (KDa):	63-68

Immunogen

Antibody Immunogen Source:	Human Plk1 (PLK) sequence peptide Cat. No.: PE-04AKY95
Antibody Immunogen Sequence:	TPN(pY)IAP(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues T214 to P220; In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Туг-217

Production

Antibody Type: Polyclonal Antibody Ig Isotype Clone Lot: Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose. Antibody Amount: 25 μg Antibody Concentration: 0.75 mg/ml Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	Troduction	
Antibody Ig Isotype Clone Lot: Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose. Antibody Amount: 25 µg Antibody Concentration: 25 µg 0.75 mg/ml Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Host Species:	Rabbit
The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose. Antibody Amount: 25 µg Antibody Concentration: 0.75 mg/ml Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Type:	Polyclonal
peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose. Antibody Amount: 25 μg Antibody Concentration: 0.75 mg/ml Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Antibody Concentration: Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Production Method:	peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification
Lot Number: 141107 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Amount:	25 μg
Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Concentration:	0.75 mg/ml
Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Lot Number:	141107
	Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
	Storage Conditions and Stability:	

AB-PK779 Plk1-pY217 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other vertebrates; Phosphosite is highly conserved in diverse species

This product is for in vitro research use only and is not intended for use in humans or animals.