

AB-PK790

Raf1-pS259 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase Raf1 (RafC)

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3



Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Raf1 proto-oncogene-encoded protein-serine kinase
Alias:	cRAF; c-RAF; c-Raf; Kinase Raf1; Murine leukemia viral (V-raf-1) oncogene homologue 1; NS5; Protein kinase raf 1; Raf proto-oncogene serine/threonine-protein kinase; RAF1; V-raf-1 murine leukemia viral oncogene 1; CCDS2612.1; ENSG00000132155; B4E1N6; Q59HC8
UniProt ID:	P04049
Sequence Predicted Mass (KDa):	75.395 (668 AA; P04049-2); 73.052 (648 AA; P04049)
Observed SDS-PAGE Mass (KDa):	70-75

Immunogen

Antibody Immunogen Source:	Human Raf1 (RafC) sequence peptide Cat. No.: PE-04ALA99
Antibody Immunogen Sequence:	RST(pS)TPN(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues R256 to N262; In the region between the C1_1 and kinase catalytic domains. This is the major in vivo phosphorylation site in Raf1.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Ser-259

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 µg
Antibody Concentration:	0.75 mg/ml
Lot Number:	141202
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerazol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK790

Raf1-pS259 Antibody



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Detection by Immunoblotting in Cell/Tissue Lysates:	Medium immunoreactivity of a target-sized protein by Western blotting in MCF7 cells with and without insulin treatment.
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins detected in MCF7 and T98G cells.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINEXUS(546-3987)