## AB-PK798 ROCK1-pY913 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase ROCK1 (ROKb)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

Tar		

Name Long:	Rho-associated protein kinase 1
Alias:	Kinase ROCK1 ;p160ROCK; Rho-associated, coiled-coil containing protein kinase 1; Rho-associated, coiled-coil forming protein kinase p160 ROCK-1; ROCK-I; ROK-beta
UniProt ID:	Q13464
Sequence Predicted Mass (KDa):	158.175 (1354 AA; Q13464)
Observed SDS-PAGE Mass (KDa):	150-170

## Immunogen

Antibody Immunogen Source:	Human ROCK1 (ROKb) sequence peptide Cat. No.: PE-04APL99	
Antibody Immunogen Sequence:	EEQ(pY)FEL(bA)C (bA) = beta-alanine	
Location in Target:	Corresponds to amino acid residues E910 to L916; Just before Rho-Binding domain. This is the major in vivo phosphorylation site in ROCK1.	
Peptide Type:	For phosphosite-specific recognition of target.	
Target Phosphosite:	Tyr-913	

## Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.	
Antibody Amount:	25 μg	
Antibody Concentration:	1 mg/ml	
Lot Number:	150305	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

## AB-PK798 ROCK1-pY913 Antibody

Applications

**Antibody Cross Reactivities:** 



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

7 (ppiloations	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended</b>	2 μg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals
<b>Antibody Positive Controls:</b>	Medium immunoreactivity with immunogen peptide on dot blots.
Detection by Immunoblotting in Cell/Tissue Lysates:	Very strong immunoreactivity of a target-sized protein by Western blotting in A431 cells, and medium target protein detection in HepG2 cells.
Overall Antibody Specificity:	Medium-very high selectivity

proteins.

In A431 cells, almost no significant cross-reactivities detected, except EGF decreases 150 KDa protein detection. In HepG2 cells, a few cross-reactive

proteins are detected, including insulin increases in 75, 50, 37 and 15 KDa

This product is for in vitro research use only and is not intended for use in humans or animals.