AB-PK803 Ros-pY2114+pY2115 Antibody

Ros-pY2114+pY2115 Antibody Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase Ros



Email: info@kinexus.ca Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Orosomucoid 1 receptor-tyrosine kinase; Proto-oncogene tyrosine-protein kinase
Alias:	ROS1 c-Ros oncogene 1, receptor tyrosine kinase; c-Ros-1; KROS; MCF3; Proto- oncogene tyrosine-protein kinase ROS; ROS1; V-ros UR2 sarcoma virus oncogene homolog 1 (avian); V-ros UR2 sarcoma virus oncogene homologue 1 (avian); MGC163394; ENSG00000169071
UniProt ID:	P08922
Sequence Predicted Mass (KDa):	263.915 (2347 AA; P08922)
Observed SDS-PAGE Mass (KDa):	200-300
Immunogen	
Antibody Immunogen Source:	Human Ros sequence peptide Cat. No.: PE-04ADN80
Antibody Immunogen Sequence:	ND(pY)(pY)RKR(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues N2112 to R2118; In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-2114+Tyr-2115
Production	Debbit
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	150203
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK803 Ros-pY2114+pY2115 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications		
Product Use:	Western blotting Antibody microarrays	
Antibody Dilution Recommended:	2 µg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals	
Overall Antibody Specificity:	Medium-High selectivity	
Antibody Cross Reactivities:	Many diverse weakly cross-reactive proteins detected in Jurkat cells.	

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)