AB-PK821 Syk-pY323 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase Syk

KINEXUS

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Spleen protein-tyrosine kinase
Alias:	DKFZp313N1010; FLJ25043; FLJ37489; Kinase Syk; KSYK; Spleen tyrosine kinase; SYK; CCDS6688.1; ENSG00000165025
UniProt ID:	P43405
Sequence Predicted Mass (KDa):	72.066 (635 AA; P43405); 69.510 (612 AA; P43405-2)
Observed SDS-PAGE Mass (KDa):	67-77

Immunogen

Antibody Immunogen Source:	Human Syk sequence peptide Cat. No.: PE-04ALE99
Antibody Immunogen Sequence:	FNP(pY)EPE(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues F320 to E326; In the region between the SH2 and kinase catalytic domains. The major site of phosphorylation in Syk.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-323

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	141107
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PK821 Syk-pY323 Antibody

Annlications

Cell/Tissue Lysates:

Overall Antibody Specificity:

Antibody Cross Reactivities:



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Very strong immunoreactivity with immunogen peptide on dot blots.
Detection by Immunoblotting in	Strong immunoreactivity of a target-sized protein by Western blotting in Jurkat

cells.

Medium-High selectivity

Almost no significant cross-reactivities detected in T98G cells and sea star oocytes, but with SDS-PAGE sample buffer extraction in Jurkat and A431 cells,

get detection of cross-reactive proteins from 120-140 KDa.

This product is for in vitro research use only and is not intended for use in humans or animals.