AB-PK834 TRIM28-pY458 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase TRIM28 (TIF1B)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

| I Or | \sim | | Jrai | מוחי |
|------|--------|---|------|------|
| 101 | ra re | - | | tein |
| | | | | |

| Name Long: | Transcription intermediary factor 1-beta | |
|--------------------------------|---|--|
| Alias: | KAP-1; KRIP-1; RING finger protein 96; RNF96; TIF1-beta; Tripartite motif-containing protein 28; TRIM28 | |
| UniProt ID: | Q13263 | |
| Sequence Predicted Mass (KDa): | 88.550 (835 AA; Q13263); 79.474 (753 AA; Q13263-2) | |
| Observed SDS-PAGE Mass (KDa): | 100-120 | |

Immunogen

| Antibody Immunogen Source: | Human TRIM28 (TIF1B) sequence peptide Cat. No.: PE-04APZ99 | |
|------------------------------|---|--|
| Antibody Immunogen Sequence: | DDP(pY)SSA(bA)C (bA) = beta-alanine | |
| Location in Target: | Corresponds to amino acid residues D455 to A461; In the region between the zf-B_box and PHD domains. One of the major in vivo sites of phosphorylation in TRIM28. | |
| Peptide Type: | For phosphosite-specific recognition of target. | |
| Target Phosphosite: | Tyr-458 | |

Production

| Antibody Host Species: | Rabbit | |
|-----------------------------------|---|--|
| Antibody Type: | Polyclonal | |
| Antibody Ig Isotype Clone Lot: | Immunoglobulin G | |
| Production Method: | The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose. | |
| Antibody Amount: | 25 μg | |
| Antibody Concentration: | 0.25 / 1.0 mg/ml | |
| Lot Number: | 150305 / 160121 | |
| Storage Buffer: | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol | |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. | |

AB-PK834 TRIM28-pY458 Antibody

Applications

Overall Antibody Specificity:

Antibody Cross Reactivities:



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

No significant cross-reactive proteins in A431, HeLa and MCF7 cells, except for

a strong immunoreactive protein in MCF7 and A431 cells.

| Applications | |
|---------------------------------------|--|
| Product Use: | Western blotting Antibody microarrays |
| Antibody Dilution Recommended: | 2 μg/ml for immunoblotting |
| Antibody Species Reactivity: | Human, mouse, rat and many other mammals |

High-very high selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.