## AB-PK839 TrkB-pY702 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-tyrosine kinase TrkB (NTRK2)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

				_					
	ar.	$\sim$	വ		r	$\neg$	t٠	$\triangle$	ın
Ta		u	U.		ш	U	ц	$\overline{}$	

Name Long:	BNDF/NT3/4/5 receptor- tyrosine kinase
Alias:	GP145-TrkB; GP145-TrkB,GP95-TrkB; NTRK2; CCDS6671.1; ENSG00000148053
UniProt ID:	Q16620
Sequence Predicted Mass (KDa):	93.826 (838 AA; Q16620-4); 91.999 (822 AA; Q16620); 81.569 (735 AA; Q16620-6); 60.994 (553 AA; Q16620-5); 59.167 (537 AA; Q16620-3); 53.051 (477 AA; Q16620-2); 35.332 (321 AA; Q16620-7)
Observed SDS-PAGE Mass (KDa):	85-95

## Immunogen

Antibody Immunogen Source:	Human TrkB (NTRK2) sequence peptide Cat. No.: PE-04AQE95
Antibody Immunogen Sequence:	RDV(pY)STD(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues R699 to D705; In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-702

## Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	150305
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## **AB-PK839** TrkB-pY702 Antibody

**Antibody Positive Controls:** 



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals; Phosphosite is highly conserved

Medium immunoreactivity with immunogen peptide on dot blots.

in mammals

This product is for in vitro research use only and is not intended for use in humans or animals.