AB-PK843 TTK-pY833+pY836 Antibody

TTK-pY833+pY836 Antibody Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase TTK



Email: info@kinexus.ca Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Dual specificity protein kinase TTK
Alias:	Cancer/testis antigen 96; CT96; ESK; Esk1; FLJ38280; Mps1; MPS1L1; Phosphotyrosine picked threonine-protein kinase; PYT; TTK protein kinase; ENSG00000112742
UniProt ID:	P33981
Sequence Predicted Mass (KDa):	97.072 (857 AA; P33981); 96.944 (856 AA; P33981-2)
Observed SDS-PAGE Mass (KDa):	90-103
Immunogen	
Antibody Immunogen Source:	Human TTK sequence peptide Cat. No.: PE-04AQH99
Antibody Immunogen Sequence:	KTL(pY)EH(pY)SGG(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues K830 to G839; In the C-terminal region after the kinase catalytic domain. These are This is the major in vivo phosphorylation sites in TTK.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-833+Tyr-836
Production Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	0.25 / 0.75 mg/ml
Lot Number:	150305 / 160121
	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Buffer:	Thosphale bunered same (TBS) priv.4, 0.05% Thimerason

AB-PK843 TTK-pY833+pY836 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Low selectivity
Antibody Cross Reactivities:	Extremely cross-reactive in MCF7 cells, but much cleaner in A431 and HeLa cells.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)