AB-PK868

MSK2-pT194 + pS196 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human protein-serine/threonine kinase MSK2 (RPS6KA4)

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3



Email: info@kinexus.ca Phone: 604-323-2547

Targ		

Name Long:	Ribosomal protein S6 kinase alpha-4
Alias:	Kinase MSK2; KS6A4; Ribosomal protein kinase B; Ribosomal protein S6 kinase, 90kDa, polypeptide 4; RPS6KA4; RSK-B; Similar to ribosomal protein S6 kinase, 90kD, polypeptide 4; Similar to ribosomal protein S6 kinase, polypeptide 4; CCDS8073.1; ENSG00000162302
UniProt ID:	O75676
Sequence Predicted Mass (KDa):	85.606 (772 AA; O75676); 84.884 (766 AA; O75676-2)
Observed SDS-PAGE Mass (KDa):	83-90

Immunogen

Antibody Immunogen Source:	Human MSK2 (RPS6KA4) sequence peptide Cat. No.: PE-04AQX85	
Antibody Immunogen Sequence:	KER(pT)F(pS)FCGC (bA) = beta-alanine	
Location in Target:	Corresponds to amino acid residues K191 to G199; In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.	
Peptide Type:	For phosphosite-specific recognition of target.	
Target Phosphosite:	Thr-194+Ser-196	

Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.	
Antibody Amount:	25 μg	
Antibody Concentration:	0.5 mg/ml	
Lot Number:	160301	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

AB-PK868 MSK2-pT194 + pS196 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.