AB-PN519 HGS-pY216 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human HGS (Hrs)

KINEXUS

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Hepatocyte growth factor-regulated tyrosine kinase substrate
Alias:	Human growth factor-regulated tyrosine kinase substrate; Protein pp110; Vps27; ZFYVE8
UniProt ID:	O14964
Sequence Predicted Mass (KDa):	86.192 (777 AA; O14964-1); 76.376 (690 AA; O14964-2)
Observed SDS-PAGE Mass (KDa):	100-110

Immunogen

Antibody Immunogen Source:	Human HGS (Hrs) sequence peptide Cat. No.: PE-04AFB95
Antibody Immunogen Sequence:	EPC(pY)EQL(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues E213 to L219; In the FYVE domain. This is the major in vivo phosphorylation site in HGS.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-216

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	0.2 / 1.0 mg/ml
Lot Number:	141017 / 160121
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-PN519 HGS-pY216 Antibody

Antibody Cross Reactivities:



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications Product Use: Western blotting | Antibody microarrays Antibody Dilution Recommended: 2 μg/ml for immunoblotting Antibody Species Reactivity: Human, mouse, rat and many other mammals

Overall Antibody Specificity:

Very high selectivity

No significant cross-reactive proteins detected in phenylarsine oxide

(PAO)+vanadate-treated HeLa cells, EGF-treated A431 cells and insulin-treated MCF7 cells, when these cells were homogenized in SDS-PAGE sample buffer.

This product is for in vitro research use only and is not intended for use in humans or animals.