## AB-PN536 SIT-pY95 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human SIT



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Name Long:	Signaling threshold-regulating transmembrane adapter 1
Alias:	gp30/40; MGC125908; MGC125909; MGC125910; RP11-331F9.5; SHP2 interacting transmembrane adaptor; SHP2-interacting transmembrane adapter protein; SHP2-interacting transmembrane adaptor protein; signaling threshold regulating transmembrane adaptor 1; Signaling threshold-regulating transmembrane adapter 1; SIT; SIT1; suppression inducing transmembrane adaptor 1; Suppression-inducing transmembrane adapter 1
UniProt ID:	Q9Y3P8
Sequence Predicted Mass (KDa):	21.126 (196 AA; Q9Y3P8)
Observed SDS-PAGE Mass (KDa):	19-21
Immunogen	
Antibody Immunogen Source:	Human SIT sequence peptide Cat. No.: PE-04AHF95
Antibody Immunogen Sequence:	NLH(pY)LQT(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues N92 to T98; In the middle of the SIT domain I the protein. One of the major in vivo phosphorylation sites in SIT.
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-95

Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	150115
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-PN536 SIT-pY95 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	Almost no significant cross-reactivities detected in HeLa cells, except phenylarsine oxide (PAO) increases 25 KDa protein; No significant cross-reactivities detected in MCF7 and T98G cells.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)