# AB-PN545 WASP-pY291 Antibody

Phosphosite-specific polyclonal antibody for monitoring the phosphorylation of human WASP



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Target Protein

Name Long:	Wiskott-Aldrich syndrome protein
Alias:	
UniProt ID:	P42768
Sequence Predicted Mass (KDa):	52.913 (502 AA; P42768-1)
Observed SDS-PAGE Mass (KDa):	54-62

#### Immunogen

Antibody Immunogen Source:	Human WASP sequence peptide Cat. No.: PE-04AFC99
Antibody Immunogen Sequence:	KLI(pY)DFI(bA)C (bA) = beta-alanine
Location in Target:	Corresponds to amino acid residues K288 to I294; In the PBD domain
Peptide Type:	For phosphosite-specific recognition of target.
Target Phosphosite:	Tyr-291

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with
Production Method:	
	0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. This antibody was also subject to negative purification over phosphotyrosine-agarose.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	150106
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-PN545 WASP-pY291 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals; Phosphosite is highly conserved in mammals
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No significant cross-reactive proteins in sea star oocytes.

This product is for in vitro research use only and is not intended for use in humans or animals.