

AB-PN571

AR-pS308 Antibody

Phosphosite-specific rabbit polyclonal antibody for AR. This phosphoserine-site antibody also cross-reacts with phosphothreonine, and glutamic acid substitutions in the phosphosite.

Website: www.kinexus.ca
Email: info@kinexus.ca
Phone: 604-323-2547



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Target Protein

Protein Name Long:	Androgen receptor
Protein Alias:	AIS; ANDR; DHTR; Dihydrotestosterone receptor; HUMARA; Kennedy disease; NR3C4; SBMA; SMAX1; testicular feminization
UniProt ID:	P10275
Protein Molecular Mass:	98,989 Da (919 AA)

Immunogen

Antibody Immunogen Source:	Synthetic phosphopeptide patterned after human AR
Antibody Immunogen Sequence:	AEY(pS)PFKGC
Antibody Immunogen Description:	Corresponds to amino acid residues A307 to G314. S310 phosphorylation inhibits transcriptional activity. This is a minor <i>in vivo</i> phosphorylation site in the androgen receptor (≥5 reports from high throughput mass spectrometry studies in PhosphoSitePlus). AR is known to be phosphorylated at this site <i>in vitro</i> by CDK1, CDK5, and CDK11B.
Antibody Target Type:	Phosphosite-specific

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Isotype:	IgG
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from these animals was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Amount:	25 µg
Antibody Concentration:	1 mg/ml
Storage Buffer:	Phosphate buffered saline pH 7.4, 0.05% Thimerasol
Storage Conditions:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Storage Stability:	>2 years

Applications

Product Use:	Western blotting Antibody microarray
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human Chimpanzee Rhesus macaque Rat Mouse

This product is for *in vitro* research use only and is not intended for use in humans or animals.

For more information on our products please visit <https://kinexus-ca.myshopify.com/> or contact us at 1-866-546-3987