

AB-PN813

CFL1-pY140 Antibody

Phosphosite-specific rabbit polyclonal antibody for Cofilin 1 (CFL1). This phosphotyrosine-site antibody is highly specific for phosphotyrosine.

Website: www.kinexus.ca
Email: info@kinexus.ca
Phone: 604-323-2547



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Target Protein

Protein Name Long:	Cofilin 1
Protein Alias:	CFL; COF1; Cofilin; Cofilin 1 (non-muscle); Cofilin, non-muscle isoform; Non-muscle isoform; P18
UniProt ID:	P23528
Protein Molecular Mass:	18,502 Da (166 AA)

Immunogen

Antibody Immunogen Source:	Synthetic phosphopeptide patterned after human CFL1
Antibody Immunogen Sequence:	CANS(pY)EEVK
Antibody Immunogen Description:	Corresponds to amino acid residues A137 to K144. The effect of Y140 phosphorylation is unclear. This is the major site of <i>in vivo</i> phosphorylation based on ≥ 1968 mass spectrometry reports recorded in PhosphoSitePlus.
Antibody Target Type:	Phosphosite-specific

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Isotype:	IgG
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from these animals was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Amount:	25 μ g
Antibody Concentration:	0.9 mg/ml
Storage Buffer:	Phosphate buffered saline pH 7.4, 0.05% Thimerasol
Storage Conditions:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at $+4^{\circ}\text{C}$ for more than 3 months. Avoid repeated freeze-thaw cycles.
Storage Stability:	>2 years

Applications

Product Use:	Western blotting Antibody microarray
Antibody Dilution Recommended:	2 μ g/ml for immunoblotting
Antibody Species Reactivity:	Human Rhesus macaque Rat Mouse

This product is for *in vitro* research use only and is not intended for use in humans or animals.

For more information on our products please visit <https://kinexus-ca.myshopify.com/> or contact us at 1-866-546-3987