

AB-PN962

L1CAM-pY1176 Antibody

Phosphosite-specific rabbit polyclonal antibody for L1CAM (CD171). This phosphotyrosine-site antibody is highly specific for phosphotyrosine.

Website: www.kinexus.ca
Email: info@kinexus.ca
Phone: 604-323-2547



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Target Protein

Protein Name Long:	Neural cell adhesion molecule L1
Protein Alias:	CAML; CAML1; CD171; CD171 antigen; HSAS; HSAS1; L1 cell adhesion molecule; MASA; MIC5; Neural cell adhesion molecule L1; S10; SPG1
UniProt ID:	P32004
Protein Molecular Mass:	140,003 Da (1,257 AA)

Immunogen

Antibody Immunogen Source:	Synthetic phosphopeptide patterned after human L1CAM
Antibody Immunogen Sequence:	TFGE(pY)RSL(βA)C
Antibody Immunogen Description:	Corresponds to amino acid residues T1172 to L1179. The effect of N-CAM L1 Y176 phosphorylation is unclear. This is a major site of <i>in vivo</i> phosphorylation in L1CAM based on ≥10 mass spectrometry reports documented in PhosphoSitePlus. L1CAM is known to be phosphorylated at this site <i>in vitro</i> by Src.
Antibody Target Type:	Phosphosite-specific

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Isotype:	IgG
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from these animals was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Amount:	25 µg
Antibody Concentration:	0.65 mg/ml
Storage Buffer:	Phosphate buffered saline pH 7.4, 0.05% Thimerasol
Storage Conditions:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Storage Stability:	>2 years

Applications

Product Use:	Western blotting Antibody microarray
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human Chimpanzee Rhesus macaque Dog Rat Mouse Chicken Frog

This product is for *in vitro* research use only and is not intended for use in humans or animals.

For more information on our products please visit <https://kinexus-ca.myshopify.com/> or contact us at 1-866-546-3987