

AB-PP546

PPP1R14B-pY29 Antibody

Phosphosite-specific rabbit polyclonal antibody for PPP1R14B (PHI-1). This phosphotyrosine-site antibody is highly specific for phosphotyrosine.

Website: www.kinexus.ca
Email: info@kinexus.ca
Phone: 604-323-2547



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Target Protein

Protein Name Long:	Protein phosphatase 1 regulatory subunit 14B
Protein Alias:	PHI-1; Phospholipase C beta 3 neighbouring gene protein; Phospholipase C beta-3 neighboring gene protein; PLCB3N; PNG; PP14B; Protein phosphatase 1, regulatory (inhibitor) subunit 14B; SOM172
UniProt ID:	Q96C90
Protein Molecular Mass:	15,911 Da (147 AA)

Immunogen

Antibody Immunogen Source:	Synthetic phosphopeptide patterned after human PPP1R14B
Antibody Immunogen Sequence:	CGPRV(pY)FQSP
Antibody Immunogen Description:	Corresponds to amino acid residues G25 to P33. The effect of Y29 phosphorylation is unclear. This is a major <i>in vivo</i> phosphorylation site in PPP1R14B based on ≥234 mass spectrometry reports recorded in PhosphoSitePlus.
Antibody Target Type:	Phosphosite-specific

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Isotype:	IgG
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from these animals was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Amount:	25 µg
Antibody Concentration:	0.3 mg/ml
Storage Buffer:	Phosphate buffered saline pH 7.4, 0.05% Thimerasol
Storage Conditions:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Storage Stability:	>2 years

Applications

Product Use:	Western blotting Antibody microarray
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human Dog Rat

This product is for *in vitro* research use only and is not intended for use in humans or animals.

For more information on our products please visit <https://kinexus-ca.myshopify.com/> or contact us at 1-866-546-3987