

# LW6273S

## Mouse Sk. Muscle Male

Mouse skeletal muscle - male - Whole lysate boiled in SDS-PAGE sample buffer



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Production

<b>Product Name Long:</b>	Mouse skeletal muscle - male - Whole lysate
<b>Production Method:</b>	Lysates were prepared from tissues that were homogenized by sonication in buffer formulated with 1% Triton X-100, 60 mM $\beta$ -glycerophosphate, pH 7.2, 20 mM MOPS, 20 mM sodium pyrophosphate, 30 mM sodium fluoride, 5 mM EDTA, 3 mM benzamidine, 2 mM EGTA, 1 mM sodium orthovanadate, 1 mM phenylmethylsulfonylfluoride, 1 mM dithiothreitol, 10 $\mu$ M leupeptin, and 5 $\mu$ M pepstatin. Lysates were prepared following sonication and 30 min ultracentrifugation at 100,000 rpm. Lysates were further diluted in SDS-PAGE sample buffer to a final concentration of 2 mg/ml.
<b>Amount:</b>	200 $\mu$ g
<b>Protein Concentration:</b>	2 mg/ml
<b>Storage Stability:</b>	1 year at -20°C

### Applications

<b>Lysate Use Description:</b>	For testing antibodies by immunoblotting.
--------------------------------	---

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-888-546-3987