

PE-01AHB95-P

KinSub1EDLYV Peptide Powder

15-mer kinase substrate peptide for assaying Lyn



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Yes-related protein-tyrosine kinase; Tyrosine-protein kinase Lyn
Name Alias:	JTK8; Kinase Lyn; KPM; LYN; V-yes-1 Yamaguchi sarcoma viral related oncogene; FLJ26625; CCDS6162.1; ENSG00000147507
UniProt ID:	P07948

Peptide Structure

Peptide Name:	KinSub1EDLYV
Peptide Origin:	KinSub1EDLYV was originally identified using a microarray with peptides that were predicted as optimal substrates for 500 human protein kinases with a proprietary algorithm developed at Kinexus with our academic partners.
Peptide Sequence Location:	Not applicable
Peptide Sequence:	FGGEEDLYVGPGGGG
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1409.5
% Peptide Purity:	> 95
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	For assaying the phosphotransferase activity of Yes-related protein-tyrosine kinase; Tyrosine-protein kinase Lyn (Lyn, UniProt ID P07948). The KinSub1EDLYV peptide demonstrated very high phosphotransferase activity with Blk, and exhibited medium specificity when assayed with over 200 other protein kinases. A listing of other kinases that show appreciable phosphotransferase activity towards this peptide are listed in Table 1.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)