

PE-01ALN95-P

KinSub8RRKSF Peptide Powder

15-mer kinase substrate peptide for assaying HIPK3



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Homeodomain-interacting protein kinase 3
Name Alias:	DJ8L15.1
UniProt ID:	Q9H422

Peptide Structure

Peptide Name:	KinSub8RRKSF
Peptide Origin:	KinSub8RRKSF was originally identified using a microarray with peptides that were predicted as optimal substrates for 500 human protein kinases with a proprietary algorithm developed at Kinexus with our academic partners.
Peptide Sequence Location:	Not applicable
Peptide Sequence:	GGRKRRKSFRRRGGG
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1730
% Peptide Purity:	> 95
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	For assaying the phosphotransferase activity of Homeodomain-interacting protein kinase 3 (HIPK3, UniProt ID Q9H422). The KinSub8RRKSF peptide demonstrated very high phosphotransferase activity with CHK2, and exhibited low specificity when assayed with over 200 other protein kinases. A listing of other kinases that show appreciable phosphotransferase activity towards this peptide are listed in Table 1.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)