

PE-01APL95-P

ESRRA (ERR1) (277-298) Peptide Powder

22-mer peptide based on ESRRA (ERR1)



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Steroid hormone receptor ERR1
Name Alias:	ERR- alpha; ERRa; ERRalpha; ESRL1; Estrogen receptor-like 1; Estrogen-related receptor alpha; Estrogen-related receptor, alpha; NR3B1
Species Origin:	Human
UniProt ID:	P11474

Peptide Structure

Peptide Name:	ESRRA (ERR1) (277-298)
Peptide Origin:	Internal sequence in target protein.
Peptide Sequence Location:	S277-R298
Peptide Sequence:	SLPLQDELAFAEDLVLDEEGAR
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Acid
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	2430.58
% Peptide Purity:	> 95
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	ULP01BAC JO_1
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)