

PE-01APS80-P

RTI Peptide Powder

3-mer peptide based on RTI



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	RTI
Species Origin:	Human

Peptide Structure

Peptide Name:	RTI
Peptide Sequence Location:	Not applicable
Peptide Sequence:	RTI
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Acid
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	388.5
% Peptide Purity:	> 80
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	USP01BMC
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)