

PE-01AQK95-P

MEK1 (333-353) Peptide Powder

22-mer immunogen peptide based on MEK1 (MKK1, MAP2K1)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	MAPK/ERK protein-serine kinase 1 (MKK1); Dual specificity mitogen-activated protein kinase kinase 1
Name Alias:	ERK activator kinase 1; Kinase MEK1; MAP kinase kinase 1; MAP2K1; MAPK,ERK kinase 1; MAPK/ERK kinase 1; MAPKK 1; MAPKK1; MKK1; PRKMK1; CCDS10216.1; Q02750; ENSG00000169032
Species Origin:	Human
UniProt ID:	Q02750

Peptide Structure

Peptide Name:	MEK1 (333-353)
Peptide Origin:	Kinase catalytic subdomain 11
Peptide Sequence Location:	E333-K353
Peptide Sequence:	EFQDFVNKCLIKNPAERADLKC
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	2566.91
Observed Peptide Mass:	2569.2
% Peptide Purity:	94.9
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	MEK1-XI
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the MEK1-XI (MEK1-5) rabbit polyclonal antibody (Cat. No.: NK099-3) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)