

PE-01ARS99-P

ErbB3 (169-183) Peptide Powder

15-mer immunogen peptide based on ErbB3 (HER3)



Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	ErbB3 (HER3) receptor-tyrosine kinase
Name Alias:	C-erbB3; ERB3; HER3; Kinase ErbB3; LCCS2; Receptor protein-tyrosine kinase erbB-3 precursor; Tyrosine kinase-type cell surface receptor HER3; V-erb-b2 erythroblastic leukemia viral oncogene 3; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; MGC88033; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; ENSG00000065361; O75810; O75811; O75812; O75813; Q9BUD7; Q9NNX2; Q9NNX3
Species Origin:	Human
UniProt ID:	P21860

Peptide Structure

Peptide Name:	ErbB3 (169-183)
Peptide Origin:	Internal sequence in target protein.
Peptide Sequence Location:	D169-C183
Peptide Sequence:	DRDAEIVVKDNGRSC
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1675.8
Observed Peptide Mass:	1674.5
% Peptide Purity:	100
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAK-18
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the ErbB3 rabbit polyclonal antibody (Cat. No.: NK231-2P) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)