

PE-01ASG99-P

A-RAF (300-314) Peptide Powder

16-mer immunogen peptide based on Raf-A (ARaf)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	A-Raf proto-oncogene serine/threonine-protein kinase
Name Alias:	ARAF; A-RAF; ARAF1; A-raf-1; Kinase A-Raf; PKS; PKS2; Proto-oncogene Pks; V-raf murine sarcoma 3611 viral oncogene; RP1-230G1.1-005; RP1-230G1_1; A6NIT1; B4DMG5; Q5H9B2; Q96II5ENSG00000078061
Species Origin:	Human
UniProt ID:	P10398

Peptide Structure

Peptide Name:	A-RAF (300-314)
Peptide Origin:	Internal sequence in target protein.
Peptide Sequence Location:	G300-K314
Peptide Sequence:	CGYYWEVPPSEVQLLK
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1911.2
Observed Peptide Mass:	1911.4
% Peptide Purity:	100
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAK-35
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the ARaf-3 rabbit polyclonal antibody (Cat. No.: NK205-4) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)