

# PE-01ASK75-P

## RAF1 (48-62) Peptide Powder

16-mer immunogen peptide based on Raf1 (RafC)



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Raf1 proto-oncogene-encoded protein-serine kinase
<b>Name Alias:</b>	cRAF; c-RAF; c-Raf; Kinase Raf1; Murine leukemia viral (V-raf-1) oncogene homologue 1; NS5; Protein kinase raf 1; Raf proto-oncogene serine/threonine-protein kinase; RAF1; V-raf-1 murine leukemia viral oncogene 1; CCDS2612.1; ENSG00000132155; B4E1N6; Q59HC8
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	P04049

### Peptide Structure

<b>Peptide Name:</b>	RAF1 (48-62)
<b>Peptide Origin:</b>	Internal sequence in target protein.
<b>Peptide Sequence Location:</b>	L48-L62
<b>Peptide Sequence:</b>	CLTDPSKTSNTIRVFL
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	None

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	1794.1
<b>Observed Peptide Mass:</b>	1794.6
<b>% Peptide Purity:</b>	76
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KMP01CAK-39
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20°C
<b>Storage Stability:</b>	Over 1 year at -20°C

### Applications

<b>Product Use:</b>	Serves as a blocking peptide for use with the Raf1-6 rabbit polyclonal antibody (Cat. No.: NK155-8) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)