

# PE-01ATX90-P

## ErbB2 (268-281) Peptide Powder

14-mer immunogen peptide based on ErbB2 (HER2, Neu)



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	ErbB2 (Neu) receptor-tyrosine kinase
<b>Name Alias:</b>	CD340; C-erbB-2; Epidermal growth factor receptor-related protein; ErbB2; HER2; HER-2; Kinase ErbB2; MLN 19; V-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian); TKR1; hCG_28177; ENSG00000141736; Q9NP09; Q9UK79
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	P04626

### Peptide Structure

<b>Peptide Name:</b>	ErbB2 (268-281)
<b>Peptide Origin:</b>	In the extracellular domain of the receptor.
<b>Peptide Sequence Location:</b>	C268-S281
<b>Peptide Sequence:</b>	CPALVITYNTDTFES
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	None

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	1559.7
<b>Observed Peptide Mass:</b>	1561.5
<b>% Peptide Purity:</b>	90
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KMP01CAK-78
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20°C
<b>Storage Stability:</b>	Over 1 year at -20°C

### Applications

<b>Product Use:</b>	Serves as a blocking peptide for use with the ErbB2-3 rabbit polyclonal antibody (Cat. No.: NK054-6) that is also available from Kinexus.
---------------------	---

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)