

PE-01AVY90-P

MEK5 (150-164) Peptide Powder

16-mer immunogen peptide based on MEK5 (MAP2K5, MKK5)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

| | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Name Long: | MAPK/ERK protein-serine kinase 5 (MKK5); Dual specificity mitogen-activated protein kinase kinase 5 |
| Name Alias: | CXA1; HsT17454; MAP kinase kinase 5; MAP2K5; MAPK,ERK kinase 5; MAPK/ERK kinase 5; MAPKK5; MKK5; PRKMK5; HsT17454; ENSG00000137764 |
| Species Origin: | Human |
| UniProt ID: | Q13163 |

Peptide Structure

| | |
|-------------------------------------|--------------------------------------|
| Peptide Name: | MEK5 (150-164) |
| Peptide Origin: | Internal sequence in target protein. |
| Peptide Sequence Location: | A150-Q164 |
| Peptide Sequence: | CAELKKILANGQMNEQ |
| Peptide N-Terminus: | Free amino |
| Peptide C-Terminus: | Amide |
| Peptide Modifications Other: | None |

Production

| | |
|-----------------------------------|-----------------------------------------------------------------------------------------|
| Peptide Production Method: | Solid-phase peptide synthesis |
| Calculated Peptide Mass: | 1789.0 |
| Observed Peptide Mass: | 1789.1 |
| % Peptide Purity: | 92 |
| Peptide Appearance: | White powder |
| Peptide Form: | Solid |
| Peptide Solubility: | Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer |
| Lot Number: | KMP01CAK-133 |
| Amount: | 1 mg |
| Storage Conditions: | Frozen at -20°C |
| Storage Stability: | Over 1 year at -20°C |

Applications

| | |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Use: | Serves as a blocking peptide for use with the MEK5-1 rabbit polyclonal antibody (Cat. No.: NK104-3) that is also available from Kinexus. |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------|

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)