

PE-01AWM99-P

AKT1 (133-147) Peptide Powder

16-mer immunogen peptide based on Akt1 (PKBa)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	RAC-alpha serine/threonine-protein kinase
Name Alias:	AKT; Akt1; AKT1 kinase; C-AKT; Kinase Akt1; PKB; PKB-alpha; PRKBA; V-akt murine thymoma viral oncogene 1; Protein kinase B; RAC; RAC-alpha serine,threonine kinase; RAC-alpha serine/threonine kinase; RAC-PK-alpha; PKBA; hCG_96740; AKT1_NEW; AKT; MGC99656; CCDS9994; B0LPE5; B2RAM5; B3KVH4; B3KXD7; B7Z5R1; Q9BV07; ENSG00000142208
Species Origin:	Human
UniProt ID:	P31749

Peptide Structure

Peptide Name:	AKT1 (133-147)
Peptide Origin:	Pre-kinase catalytic domain
Peptide Sequence Location:	E133-M147
Peptide Sequence:	CEMEVSLAKPKHRVTM
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1858.2
Observed Peptide Mass:	1857.2
% Peptide Purity:	100
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAK-148
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the Akt1-1 rabbit polyclonal antibody (Cat. No.: NK129-3) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)