

PE-01AWX85-P

AurB (319-333) Peptide Powder

16-mer immunogen peptide based on AurKB (Aurora B, AIM-1)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Aurora Kinase B (serine/threonine protein kinase 12)
Name Alias:	AIK2; AIM1; AIM-1; ARK2; AURKB; Aurora, IPL1-related kinase 2; Aurora-2; Aurora-B; STK12; AURKB; IPL1; STK12; STK5; CCDS11134.1; ENSG00000178999
Species Origin:	Human
UniProt ID:	Q96GD4

Peptide Structure

Peptide Name:	AurB (319-333)
Peptide Origin:	Kinase last alpha-chain
Peptide Sequence Location:	A319-R333
Peptide Sequence:	CAQVSAHPWVRANSRR
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1837.1
Observed Peptide Mass:	1838.8
% Peptide Purity:	85
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAK-159
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the AurKB-3 rabbit polyclonal antibody (Cat. No.: NK193-4) that is also available from Kinexus.
---------------------	---

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)