

PE-01AXL99-S

JNK1 (379-393) Peptide Solution

16-mer immunogen peptide based on JNK1 (MAPK8)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Jun N-terminus protein-serine kinase (stress-activated protein kinase (SAPK)) 1
Name Alias:	c-Jun N-terminal kinase 1; JNK; JNK1A2; JNK1-alpha-2; JNK-46;JNK21B1/2; JUN N-terminal kinase; Kinase JNK1; MAPK8; MK08; PRKM8; SAPK1; SAPK1c; Stress-activated protein kinase JNK1; CCDS7225.1; ENSG00000107643
Species Origin:	Human
UniProt ID:	P45983

Peptide Structure

Peptide Name:	JNK1 (379-393)
Peptide Sequence Location:	L379-S393
Peptide Sequence:	CLGAAVINGSQHPSSS
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1526.7
Observed Peptide Mass:	1527.2
% Peptide Purity:	100
Peptide Appearance:	Clear liquid
Peptide Form:	Solution
Peptide Solubility:	Supplied at 1 mg/ml concentration in 5% DMSO in water
Lot Number:	KMP01CAK-173
Amount:	250 µg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 6 months at -20°C and 1 mg/ml concentration

Applications

Product Use:	Serves as a blocking peptide for use with the JNK1-2 rabbit polyclonal antibody (Cat. No.: NK217-3) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)