

PE-01AZH85-P

WNK4 (1-15) Peptide Powder

16-mer immunogen peptide based on WNK4 (PRKWNK4)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Serine/threonine-protein kinase WNK4
Name Alias:	Kinase WNK4; PRKWNK4; Protein kinase with no lysine 4; Protein kinase, lysine-deficient 4; WNK lysine deficient protein kinase 4; PHA2B; CCDS11439.1; ENSG00000126562
Species Origin:	Human
UniProt ID:	Q96J92

Peptide Structure

Peptide Name:	WNK4 (1-15)
Peptide Origin:	N-terminus
Peptide Sequence Location:	M1-Q15
Peptide Sequence:	MLASPATETTVLMSQC
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1682.0
Observed Peptide Mass:	1682.7
% Peptide Purity:	88
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAK-220
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the WNK4-1 rabbit polyclonal antibody (Cat. No.: NK255-1) that is also available from Kinexus.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)