

# PE-01BFK99H-P

## Raf-B (597-603) V600E Peptide Powder

9-mer immunogen peptide based on Raf-B (BRaf)



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	B-Raf (RafB) proto-oncogene serine-threonine-protein kinase
<b>Name Alias:</b>	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; BRAF1; C-RMIL; Kinase B-Raf; P94; RAFB1; RMIL; RMIL serine/threonine-protein kinase; V-Raf murine sarcoma viral oncogene homolog B1; MGC126806; MGC138284; CCDS5863.1; ENSG00000157764
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	P15056

### Peptide Structure

<b>Peptide Name:</b>	Raf-B (597-603) V600E
<b>Peptide Origin:</b>	Internal sequence in target protein.
<b>Peptide Sequence Location:</b>	L597-R603
<b>Peptide Sequence:</b>	LATEKSR(bA)C
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	None

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	977.5
<b>Observed Peptide Mass:</b>	977.1
<b>% Peptide Purity:</b>	> 98
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KMP01CAH-1
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20°C
<b>Storage Stability:</b>	Over 1 year at -20°C

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)