

# PE-01BFM99H-P

## CTNNB1 (38-44) Peptide Powder

9-mer immunogen peptide based on CTNNB1



**KINEXUS**

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Catenin beta-1
<b>Name Alias:</b>	Beta-catenin; Catenin, beta; CATNB; CTNB1; CTNNB
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	P35222

### Peptide Structure

<b>Peptide Name:</b>	CTNNB1 (38-44)
<b>Peptide Origin:</b>	T41A mutant of beta-catenin is oncogenic in human cancers.
<b>Peptide Sequence Location:</b>	G38-P44
<b>Peptide Sequence:</b>	GATATAP(bA)C
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	None; Includes beta-alanine-cysteine at C-terminus for coupling to KLH or thio-agarose

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	761.2
<b>Observed Peptide Mass:</b>	760.8
<b>% Peptide Purity:</b>	> 98
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KMP01CAH-15
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20°C
<b>Storage Stability:</b>	Over 1 year at -20°C

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)