

# PE-01BFW99-P

## PDGFRA (839-845) Peptide Powder

9-mer immunogen peptide based on PDGFRA



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Platelet-derived growth factor receptor kinase alpha
<b>Name Alias:</b>	A1115593; Alpha platelet-derived growth factor receptor precursor; CD140a; Kinase PDGFR-alpha; MGC74795; PDGFR2; PDGFRA; PDGFR-alpha; PDGF-R-alpha; PGFRA; Platelet-derived growth factor, alpha-receptor; Rhe-PDGFR; RHEPDGFRA; CD140A; MGC74795; CCDS3495.1; ENSG00000134853
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	P16234

### Peptide Structure

<b>Peptide Name:</b>	PDGFRA (839-845)
<b>Peptide Origin:</b>	Just after kinase catalytic subdomain VII.
<b>Peptide Sequence Location:</b>	L839-H845
<b>Peptide Sequence:</b>	LARDIMH(bA)C
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	None

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	1029
<b>Observed Peptide Mass:</b>	1029.23
<b>% Peptide Purity:</b>	100
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KMP01CAL-09
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20°C
<b>Storage Stability:</b>	Over 1 year at -20°C

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)