

PE-01BHV90-P

LckSubtide Peptide Powder

15-mer kinase substrate peptide for assaying Lck



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Lymphocyte-specific protein-tyrosine kinase; Proto-oncogene tyrosine-protein kinase LCK
Name Alias:	LSK; LSK-T; p56-LCK; pp58Lck; RP4-675E8.4; LSK (T-cell specific protein-tyrosine kinase; lck tyrosine kinase; membrane associated protein tyrosine kinase; proto-oncogene LCK; put.ptk; YT16; ENSG00000182866
UniProt ID:	P06239

Peptide Structure

Peptide Name:	LckSubtide
Peptide Origin:	Developed by Kinexus based on alignment of known substrates and Kinexus Kinase Substrate Predictor v2.0 algorithm.
Peptide Sequence Location:	Not applicable
Peptide Sequence:	KKEEDDDYENLEEAC
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	None

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1900.9
Observed Peptide Mass:	1898.9
% Peptide Purity:	91.1
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP01CAR-23
Amount:	1 mg
Storage Conditions:	Frozen at -20 °C
Storage Stability:	Over 1 year at -20 °C

Applications

Product Use:	For assaying the phosphotransferase activity of Lymphocyte-specific protein-tyrosine kinase; Proto-oncogene tyrosine-protein kinase LCK (UniProt ID P06239).
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)