

PE-04AND99-P

KHS1 (171-177) pS174 Peptide Powder

9-mer immunogen and phosphatase substrate phosphopeptide based on KHS1 (MAP4K5; KHS)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Kinase homologous to SPS1/STE20 (MAP kinase kinase kinase protein-serine kinase 5 (MEKKK5))
Name Alias:	GCKR; Germinal center kinase-related; Germinal center kinase-related enzyme; KHS; Kinase KHS1; M4K5; MAP4K5; MAPK5; MAPK5; Q9Y4K4; ENSG00000012983
Species Origin:	Human
UniProt ID:	Q9Y4K4

Peptide Structure

Peptide Name:	KHS1 (171-177) pS174
Peptide Origin:	In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
Peptide Sequence Location:	K171-G177
Peptide Sequence:	KRK(pS)FIG(bA)C
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	Phosphorylated; Includes beta-alanine-cysteine at C-terminus for coupling to KLH or thio-agarose

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1088.2
Observed Peptide Mass:	1089.5
% Peptide Purity:	100
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP04CAX-41
Amount:	1 mg
Storage Conditions:	Frozen at -20°C
Storage Stability:	Over 1 year at -20°C

Applications

Product Use:	Serves as a blocking peptide for use with the KHS1-pS174 rabbit polyclonal antibody (Cat. No.: PK671) that is also available from Kinexus. This phosphopeptide may also be useful as a substrate for screening the phosphatase activity of protein phosphatases.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)