

PE-04AQF99-P

TTK (674-680) pT676 Peptide Powder

9-mer immunogen and phosphatase substrate phosphopeptide based on TTK



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Dual specificity protein kinase TTK
Name Alias:	Cancer/testis antigen 96; CT96; ESK; Esk1; FLJ38280; Mps1; MPS1L1; Phosphotyrosine picked threonine-protein kinase; PYT; TTK protein kinase; ENSG00000112742
Species Origin:	Human
UniProt ID:	P33981

Peptide Structure

Peptide Name:	TTK (674-680) pT676
Peptide Origin:	In protein kinase catalytic domain activation T-loop between subdomains VII and VIII.
Peptide Sequence Location:	D674-K680
Peptide Sequence:	DTT(pS)VVK(bA)C
Peptide N-Terminus:	Free amino
Peptide C-Terminus:	Amide
Peptide Modifications Other:	Phosphorylated; Includes beta-alanine-cysteine at C-terminus for coupling to KLH or thio-agarose

Production

Peptide Production Method:	Solid-phase peptide synthesis
Calculated Peptide Mass:	1002.0
Observed Peptide Mass:	1001.7
% Peptide Purity:	99.37
Peptide Appearance:	White powder
Peptide Form:	Solid
Peptide Solubility:	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
Lot Number:	KMP04CAX-121
Amount:	1 mg
Storage Conditions:	Frozen at -20 °C
Storage Stability:	Over 1 year at -20 °C

Applications

Product Use:	Serves as a blocking peptide for use with the TTK-pT676 rabbit polyclonal antibody (Cat. No.: PK842) that is also available from Kinexus. This phosphopeptide may also be useful as a substrate for screening the phosphatase activity of protein phosphatases.
---------------------	---

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)