

PE-04ARW90-P

AurC (185-197) pS188+pT191 Peptide Powder

12-mer immunogen and phosphatase substrate phosphopeptide based on AurKC (Aurora C, AIK3)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

| | |
|------------------------|--|
| Name Long: | Aurora Kinase C (serine/threonine-protein kinase 13) |
| Name Alias: | AIE1; AIE2; AIK3; AURKC; Aurora, Ipl1, Eg2 protein 1 and 2; Aurora, Ipl1-related kinase 3; Aurora-C; STK13; STKD |
| Species Origin: | Human |
| UniProt ID: | Q9UQB9 |

Peptide Structure

| | |
|-------------------------------------|--|
| Peptide Name: | AurC (185-197) pS188+pT191 |
| Peptide Origin: | In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII. |
| Peptide Sequence Location: | F185-R197 |
| Peptide Sequence: | FGW(pS)VHTP(pS)LPR |
| Peptide N-Terminus: | Free amino |
| Peptide C-Terminus: | Amide |
| Peptide Modifications Other: | Phosphorylated |

Production

| | |
|-----------------------------------|---|
| Peptide Production Method: | Solid-phase peptide synthesis |
| Calculated Peptide Mass: | 1716.8 |
| Observed Peptide Mass: | 1715 |
| % Peptide Purity: | 87.9 |
| Peptide Appearance: | White powder |
| Peptide Form: | Solid |
| Peptide Solubility: | Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer |
| Lot Number: | KLP04CAA-09 |
| Amount: | 1 mg |
| Storage Conditions: | Frozen at -20 °C |
| Storage Stability: | Over 1 year at -20 °C |

Applications

| | |
|---------------------|--|
| Product Use: | This phosphopeptide may be useful as a substrate for screening the phosphatase activity of protein phosphatases. |
|---------------------|--|

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINASES (546-2737)